

MOTIVATION CLIP-ON TRICYCLE

MANUFACTURER: Motivation Charitable Trust


OVERVIEW: The Clip-on Tricycle Attachment is a quick release clip-on, clip-off attachment that can be fitted to the Motivation Rough Terrain Wheelchair. It converts the wheelchair into a tricycle allowing the user to travel greater distances more effectively.

INTENDED USERS: Active adults and teens whose lifestyle includes traveling long distance or going over rough terrain.



PRODUCT CODE	WM3TRCO-01
WHEELS	20" x 1.75" (51cm x 4cm), pneumatic front wheel with 28-hole steel rim
PEDAL POSITION	The pedal position can be adjusted horizontally by approximately 5cm; vertical adjustment is fixed and is designed to suit most users
PRODUCT WEIGHT	35.3 lbs (16 kg)
BRAKES	Strong V-brake on the front wheel
CHAIN	Standard bicycle chain is included; chain tension is adjustable
GEARS	Simple single speed freewheel gear to suit most users
CERTIFICATIONS	ISO 9000
ADDITIONAL INFORMATION	Includes tricycle attachment only does not include Motivation Rough Terrain Wheelchair.

BEST FOR USER ENVIRONMENT:





9509 Vassar Avenue, Los Angeles, CA 91311, USA
TEL: +1 818.255.0100

WWW.CLASPHUB.ORG

This product information was produced by UCP Wheels for Humanity through Consolidating Logistics for Assistive Technology Supply and Provision II (CLASP II), a four-year cooperative agreement funded by the U.S. Agency for International Development under Agreement Number AID-OAA-A-17-00021, beginning August 30th, 2017.

CLASP II is a USAID-funded project implemented by UCP Wheels for Humanity (UCP Wheels) that increases access to a variety of high-quality, affordable assistive products, and promotes appropriate provision globally.